


absToolkit


—— *an open-source frequencies searcher, also a small subroutine of SIMUPKGS
distributed independently*

For reporting issue or uploading datasets, contact me via:
ykhuang@dicp.ac.cn

名称	修改日期	类型	大小
 absToolkit.exe	2020/8/17 19:55	应用程序	31,362 KB
 absToolkit_dist_accelerated.rar	2020/8/17 19:57	WinRAR 压缩文件	30,859 KB
 lib.xlsx	2020/8/17 18:27	Microsoft Excel ...	16 KB

Double click to run absToolkit.exe, remember to put lib.xlsx and absToolkit.exe in the same directory, or absToolkit will not find its dataset 😞

To get a not bad program using experience, maximize this window size. If you find trouble when want to maximize it, right click on the header of this window and, click setting to change something in it.



To choose CO absorption, press 0 and ENTER

```
will overestimate absorption energy, vdw-correction will even make problem worse. So do not blindly
trust so-called first-principle.

*Chemical environment distinguished by FWHM
full width of half max, can be used as an qualitative indicator of homogeneity of PGM_isq-CO sites
such as Ir(CO)2 on various supports, FWHM ~5cm-1, similar with organometallic complexes ~4cm-1,
indicating an nearly isolated site, see J. Phys. Chem. Lett., 2016, 7, 3854-3860.
-----
absToolkit| step 1: choose absorbate:
            [0] CO
            [1] NO
absToolkit| (to directly quit program, enter 999) >> 0
```

For instance, if you want to look up CO frequencies of moiety CO-Au/support, type-in Au and ENTER

```
-----
absToolkit| step 1: choose absorbate:
            [0] CO
            [1] NO
absToolkit| (to directly quit program, enter 999) >> 0
absToolkit| for molecule you choose, there are 37 lines of data available.
-----
absToolkit| step 2: type-in element symbol of your metal, currently supported metals: ['Rh', 'Pt', 'Au']
absToolkit| (to directly quit program, enter 999) >> Au
```

Enter frequency you measured in experiment, e.g., 2100, then ENTER

```
absToolkit| step 2: type-in element symbol of your metal, currently supported metals: ['Rh', 'Pt', 'Au']
absToolkit| (to directly quit program, enter 999) >> Au
-----
absToolkit| step 3: type-in frequency in cm-1 you want to look up.
absToolkit| >> 2100
```

Just ENTER, under most common circumstances, this function is not needed.

```
absToolkit| step 3: type-in frequency in cm-1 you want to look up.
absToolkit| >> 2100
-----
absToolkit| step 4: (press Enter directly or choose 'normal mode' if you dont know what this step means)
[1] normal mode (0)
[2] very tight (0.0001)
[3] tight (0.001)
[4] default (0.01)
[5] loose (0.1)
[6] very loose (0.3)
[7] awfully loose (0.5)
[8] Do you feel lucky? (1.0)
----->>> I just salute to Gaussian software :)
absToolkit| >> _
```

```
absToolkit| >>
-----SEARCH RESULT-----
  active metal support frequency domain mode description          reference          reference.1
11      Au      TiO2          2100 CO on anionic Au J. Phys. Chem. C 2011, 115, 22400-22408 J. Am. Chem. Soc. 2006, 128, 6341
-----
absToolkit| thank you, absorption toolkit will quit. Results above only the first two reference is listed.
If you need more reference to cite, search the number emerges at head of line, it is the line number
of datafile lib.xlsx.

If you want to have contribution on this program,
-> send your customized lib.xlsx to me ykhuang@dicp.ac.cn
-> make comments and report issues at https://github.com/kirk0830/SIMUPKGS

absToolkit| press ENTER to quit. >> _
```

Finally we obtain relating data (if exists), “mode description” is a short description of vibration, also there are two references provided explicitly. For more information, please read instructions of absToolkit.

